

KONICA CR TECHNIQUE CHART 66% More mAs

Part	View	Small		Medium		Large	
		kV	mAs	kV	mAs	kV	mAs
Abdomen	AP (Grid) Most Dose	80	11 - 13	80	24 - 28	80	37 - 50
Ankle	AP	65	2.4	65	2.7	65	3.1
Ankle	Obl	65	2.1	65	2.5	65	2.9
Ankle	Lat	65	2	5	2.3	5	2.7
Chest -Adult	AP (non grid - 72")	85	1.9 - 2.4	85	2.9 - 3.3	90	4.3 - 5.3
Chest -Adult	Lat (non grid - 72")	90	3.3 - 4.8	90	6 - 8.4	90	10.7 - 13.3
Chest - Baby	PA (72")	80	2.4	80	3.6	80	4.4
Chest - Baby	Lat (72")	80	5.7	80	7.1	80	8.9
Chest - Newborn	AP (40")	70	2.1	72	1.6	74	1.8
Chest - Newborn	Lat (40")	74	2.8	76	3	78	3.4
C-Spine	AP (Bucky - 72")	80	18	80	25	80	2-Feb
C-Spine	AP (Bucky - 40")	80	6.3	80	8	80	10.7
C-Spine	Odontoid (72")	80	23.3	80	32	80	43
C-Spine	Odontoid (40")	80	7.4	80	10.7	80	14
C-Spine	Lat (Bucky - 72")	80	20 - 25	80	26 - 32	80	36 - 43
C-Spine	Swimmers (40")	90	36 - 57	95	57 - 69	100	80 - 100
C-Spine	Trauma Obl. (non grid)	70	14.7	77	20	77	25
C-Spine	AP (40")	70	4.4	70	5.3	73	6.7
C-Spine	Lat (72")	70	14.7	70	17.7	70	21.3
Elbow	AP	65	2.9	65	3.2	65	3.6
Elbow	Obl	65	3.2	65	3.6	65	3.9
Elbow	Lat	65	2.8	65	3.2	65	3.6
Femur - Distal	Lateral (non grid)	70	3.9	70	4.4	70	5.7
Finger	All Views	63	1.5	63	1.7	63	2.1
Foot	AP	65	2.1	65	2.8	65	3.6
Foot	Obl	65	2.5	65	3.2	65	3.9
Foot	Lat	65	3.2	65	3.9	65	4.9
Forearm	AP	65	2.1	65	2.5	65	2.8
Forearm	Lat	65	2.3	65	2.7	65	4.9
Hand	PA	65	1.7	65	2.1	65	2.5
Hand	Obl	65	2.1	65	2.5	65	2.8
Hand	Lat	65	2.5	65	3.2	65	3.9



Digital Radiography Solutions "Maximizing Image Quality, Minimizing Patient Dose"

Disclaimer: Digital Radiography Solutions (dRs) cannot be held responsible for any images taken or passed at your facility. Your radiologist or physician always has the final say if the image is of diagnostic value.