

KONICA CR TECHNIQUE CHART Least mAs

Part	View	Small		Medium		Large	
		kV	mAs	kV	mAs	kV	mAs
Abdomen-(LgM 2.1)	AP (Grid)	80	12 -16	80	24 - 28	80	35 - 45
Ankle	AP	65	1.3	65	1.5	65	1.7
Ankle	Obl	65	1.2	65	1.4	65	1.6
Ankle	Lat	65	1.1	5	1.3	5	1.5
Chest -Adult	AP (non grid - 72")	85	1 - 1.2	85	1.6 - 2	90	2.5 - 3.2
Chest -Adult	Lat (non grid - 72")	90	2.2 - 2.8	90	3.6 - 4.5	90	6.3 - 8
Chest - Baby	PA (72")	80	1.6	80	2	80	2.5
Chest - Baby	Lat (72")	80	3.2	80	4	80	5
Chest - Newborn	AP (40")	70	1.2	72	1.2	74	1.2
Chest - Newborn	Lat (40")	74	2	76	2	78	2
C-Spine	AP (Bucky - 72")	80	10	80	14	80	18 to 20
C-Spine	AP (Bucky - 40")	80	3.5	80	4.5	80	6
C-Spine	Odontoid (72")	80	12 to 14	80	17-18	80	24
C-Spine	Odontoid (40")	80	4.2	80	5.5	80	7
C-Spine	Lat (Bucky - 72")	80	11 to 14	80	15-18	80	20 - 24
C-Spine	Swimmers (40")	90	20 - 32	90	32 - 40	90	40 - 50
C-Spine	Trauma Obl. (non grid)	70	8	70	11	70	14
C-Spine	Lat (non grid - 72")	70	8	70	10	73	12.0 - 13.0
C-Spine	AP (non grid - 40")	70	2.5	70	3	70	3.5 - 4
Elbow	AP	65	1.6	65	1.8	65	2
Elbow	Obl	65	1.6	65	1.8	65	2
Elbow	Lat	65	1.6	65	1.8	65	2
Femur - Distal	Lateral (non grid)	70	1.8	70	3.2	70	4
Finger	All Views	63	0.6	63	0.8	63	1
Foot	AP	65	1.5	65	1.7	65	2
Foot	Obl	65	1.7	65	2	65	2.5
Foot	Lat	65	1.8	65	2.4	65	2.8
Forearm	AP	65	1.5	65	1.6	65	1.7
Forearm	Lat	65	1.6	65	1.7	65	1.8
Hand	PA	65	0.7	65	0.9	65	1.1
Hand	Obl	65	0.8	65	1.1	65	1.3
Hand	Lat	65	1.5	65	1.8	65	2.2



Digital Radiography Solutions "Specializing in Educational Conferences and Seminars"

Disclaimer - Digital Radiography Solutions (DRS) cannot be held responsible for any images taken or passed in your facility. Your radiologist always has the final say if an image is diagnostic.