

## KONICA CR TECHNIQUE CHART      Least mAs

Part	View	Small		Medium		Large	
		kV	mAs	kV	mAs	kV	mAs
Hip	AP ( non grid )	70	3.2	70	4	70	5
Hip	X-Table Lat (Grid)	90	16 - 25	90	30 - 40	90	50 - 60
Humerus	AP ( non grid )	65	1.25	65	2	65	2.5
Knee	AP (Bucky)	75	2.7	75	3.6	75	4
Knee	Obl (Bucky)	75	2.4	75	3.2	75	3.6
Knee	Lat (Bucky)	75	2.4	75	3.2	75	3.6
Knee	Sunrise (non grid)	65	3	65	3.8	65	4.8
Knee	Non-Bucky	65	2.7	65	3	65	3.3
L-Spine	AP	80	7 to 11	80	14 to 18	80	22 - 26
L-Spine	X-Table Lat (Grid)	90	55 - 65	90	85 - 100	90	160 - 200
Mandible	Obl (non grid - 40")	65	4.5	65	6	65	7.5
Pelvis	AP (Grid)	80	7.3	80	14.5	80	23
Ribs	Upper (72")	75	8 to 12	75	14 to 20	75	25 - 32
Ribs	Lower (40")	80	5.0 - 8	80	10 - 12.5	80	16 - 20
Ribs	Obl Upper (72")	75	7 to 9	75	14 - 22	75	23 - 28
Shoulder	AP (non grid)	65	2.4	65	3	65	3.2
Shoulder	Mercedes (Bucky)	75	7.5	75	ten	75	fifteen
Shoulder	Axillary (non grid)	65	3.7	65	4.5	65	6
Sinus	Caldwell	80	6	80	7.5	80	9
Sinus	Waters	80	7.5	80	9	80	10
Sinus	Lateral	80	3	80	3.7	80	4.8
Skull	AP	80	9	80	12	80	13.5
Skull	Lat (Grid)	80	3.7	80	4.8	80	5.5
Tib-Fib	AP (tt)	65	2.4	65	3	65	4
Tib-Fib	Lat (tt)	65	2	65	2.4	65	3
Toe	All Views	60	1.7	60	2.2	60	2.8
T-Spine	AP	80	6.5 - 9	80	15 - 18	80	28
T-Spine	Lat (2 sec)	85	11 to 16	85	25 - 28	85	45 - 50
Wrist	PA	60	1.3	60	1.5	60	1.7
Wrist	Obl	60	1.5	60	1.7	60	1.9
Wrist	Lat	65	1.5	65	1.7	65	1.9
Zygomatic Arch	SMV view (non grid-30")	65	7.5	65	9	65	11



Digital Radiography Solutions "Specializing in Educational Conferences and Seminars"

**Disclaimer** - Digital Radiography Solutions (DRS) cannot be held responsible for any images taken or passed in your facility. Your radiologist always has the final say if an image is diagnostic.