

## UNIVERSAL CR TECHNIQUE CHART 66% More mAs

Part	View	Small		Medium		Large	
		kV	mAs	kV	mAs	kV	mAs
Hip	AP ( non grid )	77	3.6	77	4.4	77	5.7 - 6.9
Hip	X-Table Lat (Grid)	90	28.4 - 44	90	53 - 69	90	90 - 107
Humerus	AP	70	2.1	70	3.6	70	4.4
Knee	AP (Bucky)	81	3.2	85	4.4	85	5.7
Knee	Obl (Bucky)	81	3.2	85	4.4	85	5.7
Knee	Lat (Bucky)	81	2.8	85	3.6	85	4.4
Knee	Sunrise ( non grid )	70	3.6	70	4.4	70	5.7
Knee	Non-Bucky	70	3.2	70	3.9	70	4.4
L-Spine	AP	90	7.1 - 11.2	90	14.1 - 17.7	90	22.3 - 28.4
L-Spine	X-Table Lat (Grid)	95	71 - 89	95	107 - 143	95	177 - 284
Mandible	Obl (non grid - 40")	70	5.2	70	6.9	70	8.9
Pelvis	AP (Grid)	85	8.9	85	16	85	24.9
Ribs	Upper (72")	80	6.9 - 11.2	80	12.4 - 17.7	80	22.3 - 28.3
Ribs	Lower (40")	85	8.9 - 14.1	85	17.7 - 24.9	85	28.4 - 35.5
Ribs	Obl (72")	80	8.9 - 17.7	80	21.3 - 32	80	36 - 49.3
Shoulder	AP (non grid)	77	2.8	77	3.9	77	4.9
Shoulder	Y View	81	8.9	81	12..8	81	17.7
Shoulder	Axillary (non grid)	70	4.4	77	5.7	77	6.9
Sinus	Caldwell	85	6.9	85	8.9	85	11.2
Sinus	Waters	85	8.9	85	11.2	85	12.4
Sinus	Lateral	85	3.6	85	4.4	85	5.7
Skull	AP	85	8.9	85	11.2	85	14.3
Skull	Lat (Grid)	85	4.4	85	5.7	85	7.1
Tib-Fib	AP	77	2.8	77	3.3	77	3.9
Tib-Fib	Lat	77	2.5	77	2.9	77	3.6
Toe	All Views	63	2.1	63	2.7	63	3.6
T-Spine	AP	90	6.4 - 8.9	90	14.3 - 17.7	90	24.9 - 28.4
T-Spine	Lat (2 sec)	90	14.3 - 22.3	90	28.4 - 36	90	44 - 64
Wrist	PA	66	1.5	66	1.7	66	2.1
Wrist	Obl	66	1.6	66	2	66	2.5
Wrist	Lat	70	8.9	70	10.7	70	12.9
Zygomatic Arch	SMV view (non grid - 30")	70		70		70	



Digital Radiography Solutions "Specializing in Educational Conferences and Seminars"

**Disclaimer** - Digital Radiography Solutions (DRS) cannot be held responsible for any images taken or passed in your facility. Your radiologist always has the final say if an image is diagnostic.