

Difference Between All Technique Charts (in Speed and mAs)								
DR Cesium Least mAs	DR Cesium 33% More	DR Cesium 66% More	DR Cesium 100% More	DR Gado- linium	CR Least mAs	CR 33% More	CR 66% More	CR 100% More
600 Speed	500 Speed	400 Speed	300 Speed	450 Speed	450 Speed	375 Speed	300 Speed	225 Speed
10 mAs	13.3 mAs	16.6 mAs	20 mAs	15 mAs	15 mAs	20 mAs	25 mAs	30 mAs



Digital Radiography Solutions “Maximizing Image Quality, Minimizing Patient Dose”

Disclaimer: Digital Radiography Solutions (dRs) cannot be held responsible for any images taken or passed at your facility. Your radiologist or physician always has the final say if the image is of diagnostic value.