

CR EI Comparisons (Agfa/Fuji/Carestream)

Agfa - LgM	Fuji - S	Carestream - EI
1.80	400.0	1700
1.81	393.3	1710
1.82	386.7	1720
1.83	380.0	1730
1.84	373.3	1740
1.85	366.7	1750
1.86	360.0	1760
1.87	353.3	1770
1.88	346.7	1780
1.89	340.0	1790
1.90	333.3	1800
1.91	326.7	1810
1.92	320.0	1820
1.93	313.3	1830
1.94	306.7	1840
1.95	300.0	1850
1.96	293.3	1860
1.97	286.7	1870
1.98	280.0	1880
1.99	273.3	1890
2.00	266.7	1900
2.01	260.0	1910
2.02	253.3	1920
2.03	246.7	1930
2.04	240.0	1940
2.05	233.3	1950
2.06	226.7	1960
2.07	220.0	1970
2.08	213.3	1980
2.09	206.7	1990

Agfa - LgM	Fuji - S	Carestream - EI
2.10	200.0	2000
2.11	196.7	2010
2.12	193.3	2020
2.13	190.0	2030
2.14	186.7	2040
2.15	183.3	2050
2.16	180.0	2060
2.17	176.7	2070
2.18	173.3	2080
2.19	170.0	2090
2.20	166.7	2100
2.21	163.3	2110
2.22	160.0	2120
2.23	156.7	2130
2.24	153.3	2140
2.25	150.0	2150
2.26	146.7	2160
2.27	143.3	2170
2.28	140.0	2180
2.29	136.7	2190
2.30	133.3	2200
2.31	130.0	2210
2.32	126.7	2220
2.33	123.3	2230
2.34	120.0	2240
2.35	116.7	2250
2.36	113.3	2260
2.37	110.0	2270
2.38	106.7	2280
2.39	103.3	2290
2.40	100.0	2300



Digital Radiography Solutions "Specializing in Educational Conferences and Seminars"

Disclaimer - Digital Radiography Solutions (DRS) cannot be held responsible for any images taken or passed in your facility. Your radiologist always has the final say if an image is diagnostic.