

## CR/DR EI Comparisons (Konica/Shimadzu/SwissRay)

S (Konica)	EXI (Shimadzu)	DI (Swissray)
400	200	25.0
395	205	25.5
390	210	26.0
385	215	26.5
380	220	27.0
375	225	27.5
370	230	28.0
365	235	28.5
360	240	29.0
355	245	29.5
350	250	30.0
345	255	30.5
340	260	31.0
335	265	31.5
330	270	32.0
325	275	32.5
320	280	33.0
315	285	33.5
310	290	34.0
305	295	34.5
300	300	35.0
295	305	35.5
290	310	36.0
285	315	36.5
280	320	37.0

S (Konica)	EXI (Shimadzu)	DI (Swissray)
275	325	37.5
270	330	38.0
265	335	38.5
260	340	39.0
255	345	39.5
250	350	40.0
245	355	40.5
240	360	41.0
235	365	41.5
230	370	42.0
225	375	42.5
220	380	43.0
215	385	43.5
210	390	44.0
205	395	44.5
200	400	45.0
190	440	47.5
180	480	50.0
170	520	52.5
160	560	55.0
150	600	57.5
140	640	60.0
130	680	62.5
120	720	65.0
110	760	67.5
100	800	70.0



Digital Radiography Solutions *“Specializing in Educational Conferences and Seminars”*

**Disclaimer** - Digital Radiography Solutions (DRS) cannot be held responsible for any images taken or passed in your facility. Your radiologist always has the final say if an image is diagnostic.